

## Installing the USB Driver for the Reference Light Meter


The first time the Reference Light Meter is plugged into a USB port on your computer, Windows will need to search for the appropriate USB driver. These are located in the Chroma5 USB Driver folder.

If you have plugged in the Reference Light Meter for the first time, you will need to install the Windows driver. The following window appears.




## Installing the USB Driver for the Reference Light Meter

- 1 Select the NO, NOT THIS TIME radio button and then click NEXT.



- 2 Select the INSTALL FROM A LIST OR SPECIFIC LOCATION (ADVANCED) radio button and click NEXT.
- 3 Check the INCLUDE THIS LOCATION IN THE SEARCH box.



- 4 Make sure the folder location below it is pointing to the Chroma5 USB Driver folder. Click NEXT to continue.



- 5 For a few seconds, you will see a dialog window as the driver is installed. When it is complete, the last window appears. Click FINISH to close the window.